

Spectroscopic lamps



Области применения

- Laboratory & Analysis
- Spectroscopy


CBP_09.00 Product features and benefits, all [calc.]

- High quality monochromatic light sources
- Intense and stable irradiance at specified line spectra
- Suitable for AC operation




(Hg) CONTAINS MERCURY

Product family datasheet



SPECTRAL LAMPS



SPECTRAL LAMPS

Product family datasheet

Technical data

Product description	General Product Information	Electrical Data	
		Lamp voltage	Lamp current
Cd/10	Cadmium	15 V	1. A
Cs/10	Caesium	10 V	1. A
He/10	Helium	60 V	1. A
Na/10	Sodium	15 V	1. A
Ne/10	Neon	30 V	1. A
Hg 100	Mercury	45 V	1. A
HgCd/10	Mercury/Cadmium	30 V	1. A
Tl/10	Thallium	15 V	1. A
Zn/10	Zinc	15 V	1. A

Product family datasheet

Application advice

For more detailed application information and graphics please see product datasheet.

Disclaimer

Subject to change without notice. Errors and omission excepted. Always make sure to use the most recent release.